

The Analytik Jena qTOWER³ is Now Available for Demo!

The qTOWER³ product family guarantees well-founded real-time results. It benefits from 30-years of experience, adapting peerless temperature control precision in the sample block whether it be up to 96 or up to 384 samples.

The patented high-performance optics guarantee the outstanding homogeneous excitation and illumination of all individual samples within 6 sec per block. The qTOWER³ product family achieves unique flexibility with its proven filter module equipment, which can be freely configured and expanded at any time and enables up to six channel multiplexing. Find out how you can benefit with Analytik Jena qTOWER³ systems today.

- High-end real-time PCR thermal cycler with quick ramp rates and true gradient function
- Multiplex analysis: freedom to equip up to 6 excitation and emission filters (field installable)
- Large selection of modules to choose from
- Stand Alone System controlled by in-built touch screen tablet – analysis on board or off-line
- Industry Leading 10-year warranty on optical components (2-year on system – extendable)

Analytik Jena qTOWER³ G 96 Touch qPCR system

- 10 inch colour touch tablet for setup and analysis
- qPCRsoft package for convenient operation and control
- User Extendable Filter Module System for maximum flexibility
- Patented Fiber Optic System for ideal real-time PCR signals
- High-quality silver block for maximum thermal conductivity


Vision, amazement, simplicity:

Get in touch with Bio-Strategy to trial the qTOWER³ real-time PCR thermal cycler and complete your workflow.

Description	Code
qTOWER ³ G touch (stand-alone device, no need for PC) with gradient and colour module 1 - FAMTM, Sybr [®] Green, Alexa488	ANJE844-00556-2
Available Colour Modules to trial!	
Colour module 2 – JOETM, HEXTM, VIC [®] , YakimaYellowTM	ANJE844-00521-0
Colour module 3 – TAMRATM, DFOTM, Alexa546 [®] , NEDTM	ANJE844-00522-0
Colour module 4 – ROXTM, TexasRed [®] , Cy3,5 [®]	ANJE844-00523-0
Colour module 5 – Cy5 [®] , Alexa633 [®] , Quasar670TM	ANJE844-00524-0



Ask us to trial the qTOWER³ Touch system today!



Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Rev. April 2021

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.