

Epredia SlideMate Pro Available from Bio-Strategy!

Dual-hoppers. Singular solution.

- **Efficient performance: dual-hopper design with intelligent slide selection**
- **High throughput: 4-5 seconds per slide**
- **Simple integration: Windows based driver and Gateway functions**

Introducing the SlideMate Pro thermal transfer microscope slide printer from Epredia, a smart solution with flexible configuration options to match your laboratory's needs.

Avoid the annoyance and wasted time caused by constantly loading and unloading your laboratory's slide printer to switch between multiple slide types. The SlideMate Pro printer features a dual-hopper loading system with intelligent slide selection.

Load adhesion slides in one hopper, non-adhesion slides in the other – the printer can be programmed to automatically read the requested slide type from the incoming print image and pull from the appropriate hopper.

Have a need for one-off or specialty slides? Simply load them through the manual feed as required – no need to remove the hoppers.

Want to maximise throughput? Load matching slides in both hoppers and enjoy the benefits of quick printing times.

Setup is simple. The new SlideMate Pro software uses a Windows® based driver that easily connects to other computers and laboratory information systems. SlideMate Pro enables sample-tracking solutions including user log-in information, time and date stamp on each slide, and unique transaction IDs.

Introductory Price!
Contact Us For More Details!



Product Specifications

Dimensions (DxWXH)	330 x 163 x 406 mm 13 x 6.4 x 16 in
Power Requirements	Universal 100-240 VAC 50-60 Hz, 1.2A Max.
Slide Type	25 x 75 x 1 mm, 1 x 3 in
Print Speed	4-5 seconds per slide
Media Type	Thermal transfer ribbon
Media Capacity	1 roll, ~12,000 slides per roll
Media Colour / Resolution	Black / 300 dpi
Connectivity	Ethernet, USB-A (2), USB-B (1)

Description

Description	Code
SlideMate Pro Printer	EPRNA83900015
Spare Hopper	EPRNA83910011
Thermal Transfer Ribbon	EPRNA83910001
Head Cleaning Kit	EPRNA83910008

POWERED BY

SHANDON

MICROM

MENZEL-GLÄSER

Richard-Allan Scientific

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2021

Ref. Code:
EPRN-MATE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.