

Hach DR890 Colorimeter

Standard features of the DR890 Portable Colorimeter Series with 90 methods and automatic wavelength selection.

Please, discuss your application with us and we will help you to select the appropriate consumables.

- Pre-programmed: ready-to-use out of the box
- Compact: true hand-held analysis, for use anywhere
- Rugged: dustproof, waterproof, shockproof, worry-free (designed to meet IP67 standards)
- Powerful: full datalogging and RS232 interface capabilities
- Simple: menu-driven, step-by-step analysis
- Affordable: choose the instrument you need or select a complete portable laboratory batteries-/accu operation possible, Real-Time-Clock, integrated timer, data memory with 100 measuring values and automatic AutoOff-Function, wavelengths: 420nm, 520nm, 560nm and 610nm
- Easy-to-use software, pre-programmed methods, and prepared reagents
- Designed for use in the field
- Rugged construction
- Reliable results without a mains connection
- Configured for immediate use

Each DR/800 Series Colorimeter includes two 1-in. round sample cells marked at 10, 20 and 25ml, two 1cm sample cells, 16mm COD/TNT adapter, batteries and illustrated instrument and procedure manuals.

On Sale!
(While Stocks Last)



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
DR890 Colorimeter , 1-Month Warranty	HACH48470-00	N/A	\$1,457
HQ11D Portable pH Meter , 4 AA batteries, PHC201 Laboratory Gel pH Electrode, 1m Cabl, Intellical Probe Stand, HQD Meter Stand, pH Buffer Solutions, 4.01/7.00/10.01, 500ml Bottles, pH Storage Solution, 500ml Bottle, Quick-Start Guide, Test Certificate and User Manuals	HACH8507530	\$2,027	\$1,946
Conductivity Starter Kit : HQ14d Conductivity Meter, Standard Conductivity Probe with 1m Cable, HQd Meter Stand, Conductivity Standard (1000µS/cm, 100ml), and Probe Stand	HACH8506100	\$2,306	\$2,213
Complete Water Quality Lab Kit : DO/pH/Conductivity, HQ40d Meter with A/C Power and USB Adapters, Refillable pH Probe with 1m Cable, Conductivity Probe with 1m Cable, Standard LDO Probe with 1m, Cable, HQd Meter Stand, pH Buffer Kit (500ml), pH Storage Solution (500ml), Conductivity Standard, (1000µs/cm, 100ml), Disposable BOD Bottle and Stopper for Calibration, and Probe Stand	HACH8506000	\$5,505	\$5,285
DRB200 Reactor for Digestion of up to 30 COD/TNT Vials	HACHLTV082.52.44001	\$2,721	\$2,612
Digital Titrator	HACH1690001	\$340	\$326

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2021 or
While Stocks Last

Ref. Code:
HACH-COLOR

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.