

## Heidolph Mr Silver Package: Magnetic Stirrer with Precise Temperature Control

### External temperature sensor for precise control.

- Speed range: 100 to 1,400 rpm
- Mixes up to 20L
- Max. temperature of 300°C
- Chemically resistant, scratch-proof
- Durable construction: 10-year lifespan
- Precise temperature control
- 3-year warranty

Description	Code	Price	Sale
Heidolph MR Silver 1 Package Includes Magnetic Stirrer MR Hei-Tec, Temperature Sensor Pt 1000 (V4A), and Clamping System	HEID505-30080-00	\$1,457	<b>\$1,385</b>

The standard RS 232 interface of the MR Hei-Connect and MR Hei-End enables the **documentation and controlling via PC**. A free software is available on our website

The extended heating capacity of **800 W reduces heat-up times by 35 %** compared to other units at 600 W

The chemically resistant Kera-Disk® hot plate allows for **immediate heat transfer** resulting in quick heat-up times

**Hermetically-sealed housing protects** all mechanical and electronic components **from aggressive environments**

In case the heating function fails, stirring will not be discontinued to **prevent bumping**

**Save Today!**  
(While Stocks Last)

The temperature sensor consists of **two independent safety circuits** which switches off heating in case of any overtemperature situation

In case of a short-circuit, a damage or removal of the temperature sensor from the media, the **unit powers off completely**

Damage to the stirrer is categorically ruled out even if exposed to highest temperatures – all models come with a **fire-resistant aluminum die-cast housing**

A separate on/off button for heating **prevents unintentional heat-up** – the button is illuminated for visual control. In addition, the residual heat indicator **prevents from potential burning hazards**



Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453  
E: sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 June 2021 or  
While Stocks Last!

Ref. Code:  
HEID-HEITEC

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.