

Kimble Blood and Haematocrit Tubes

Kimble tubes for blood typing and Micro-Haematocrit tubes for checking the red blood cell percentage in blood.

Blood Tubes:

- Kimble borosilicate blood typing tubes feature a permanent coloured ceramic label, pre-labelled with blood type

Haematocrit Tubes:

- Kimble micro-haematocrit tubes supplied in reclosable snap-lock vials to maintain cleanliness
- Ends may be sealed by flame or sealing wax/clay
- Non-heparin or heparinised to prevent blood clotting
- Ammonium heparin does not interfere with sodium identification
- Manufactured from soda-lime glass



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Tube Blood Anti A Blue 3ml 10x75mm Borosilicate (1,000)	KIML60A10BZ2	\$162.90	\$122.00
Tube Blood B Blue 3ml 10x75mm Borosilicate (1,000)	KIML60A10BZ3	\$162.90	\$122.00
Tube Blood Anti B Yellow 3ml 10x75mm Borosilicate (1,000)	KIML60A10BZ4	\$162.90	\$122.00
Tube Blood Is Green 3ml 10x75mm Borosilicate (1,000)	KIML60A10BZ5	\$162.90	\$122.00
Tube Blood Iis White 3ml 10x75mm Borosilicate (1,000)	KIML60A10BZ6	\$162.90	\$122.00
Tube Mark-M White 3ml 10x75mm Borosilicate (1,000)	KIML60A10BZW	\$174.80	\$131.00
Tube Blood Anti A Blue 5ml 12x75mm Borosilicate (1,000)	KIML60B12BZ2	\$179.00	\$134.00
Tube Blood Anti B Yellow 5ml 12x75mm Borosilicate (1,000)	KIML60B12BZ4	\$181.60	\$136.00
Tube Mark-M White 5ml 12x75mm Borosilicate (1,000)	KIML60B12BZW	\$178.20	\$133.00
Tube Haematocrit Mini Capillary Amm Heparin Black (100/Vial)	KIML40C505	\$15.00	\$11.00
Tube Haematocrit Micro Capillary Non-Heparin Blue (200/Vial)	KIML41A2502	\$17.30	\$13.00
Tube Haematocrit Micro Capillary Amm Heparin Red (200/Vial)	KIML41B2501	\$15.00	\$11.00

For more available tubes from DWK, contact Bio-Strategy today!

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2021

Ref. Code:
KIML-BLOOD

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.