

Miele PG8583 Washer/Disinfector

Free installation and commissioning*, free detergent and neutraliser (total savings: \$942).

Package includes A102 upper basket with 2 x AK12 inserts, A150 lower module basket with 2 x A301/5 modules (3 x 6 = 18 nozzles/module), and waterproofing kit.

This underbench /standalone model combines a powerful pump with integrated heating elements for safe washing cycles up to 93°C . The EcoDry passive drying automatically opens the door at the end of the run when temperature has dropped below 70°C allowing the load to dry and cool down faster.

- 15 programs including new “Mini” program for lightly soiled glassware
- Configured for washing wide neck and narrow neck glassware
- Dispensing pumps for liquid detergent and liquid neutraliser
- Compact: 600mm x 600mm x 850mm (820mm with lid removed)
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Many other options and accessories to select from
- Spray arm and spray pressure monitoring

Reduce the risk of injury through manual washing and improve safety in your wash area!

Package Deal!



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Miele PG8583 Washer/Disinfector Complete Package 1	MIEL PG8583	\$13,675	\$12,733

**(city metro area's only within 100km of GPO).*

Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2021

Ref. Code:
MIEL-PG8583

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.