

Priorclave Q63 Series Autoclaves

Large capacity front loading cylindrical chamber autoclaves.

Impressive load capacities ideally suited for sterilising larger loads and bulky items.

- Available in capacities of 320 or 400 litres with 630mmØ chamber
- Quickseal single-action door mechanism fitted with pressure and thermal locks
- Grade 316 stainless steel chamber carries a 10-year warranty
- Many options including pre-cycle vacuum, steam heating, trolley, loading cassette, etc
- Supplied with automatic water-feed, 2 shelves, media warming, timed free steaming, fan assisted cooling and USB port
- Low water level sensor
- Tactrol microprocessor control with simple fully variable setting of process time and temperature with graphic indication of cycle status
- Single flexible parameter program with option for 5 or 10 program capacity
- Operating range: up to 140°C; 2.4 bar



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
320 Litre Priorclave, Electrically Heated with Elements in Chamber	PRIOPS/Q63/EH320	\$43,642	\$39,278
400 Litre Priorclave, Electrically Heated with Elements in Chamber	PRIOPS/Q63/EH400	\$45,993	\$41,395

Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



VALID UNTIL:
30 June 2021

Ref. Code:
PIRO-Q63

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.